


# Gas separation unit

## Overview & assembly manual



## Notes

This electrolyser gas separation unit is intended for research purposes only and can be used for many different purposes.


There is no warranty on performance, corrosion or lifetime on the items. It is purely for research purposes.

**IMPORTANT NOTE ABOUT GOOD SAFETY PRACTICE WHEN OPERATING ELECTROLYSERS THAT GENERATE H<sub>2</sub> AND O<sub>2</sub>:** No matter if two separation units or bottles are used for separating out H<sub>2</sub> and O<sub>2</sub> on each side, there will always be a transfer of gases from one side to the other (either through solubility in water if the water from both sides is mixed at some point or by cross-over through the membrane/separator). This has the consequence that the head space in bottles with time can become a mix of both H<sub>2</sub> and O<sub>2</sub> that can reach explosive limits. For this reason it is strongly recommend to:

- Keep the head-space in bottles as low as possible - if an explosion occurs the gas volume is as small as possible
- Ensure proper ventilation and suction from gas separation units/bottles
- Keep setup in enclosures to prevent personal injuries if either a tube burst or an explosion happens - **HOWEVER IT IS STRICTLY NEEDED TO HAVE VENTILATION IN THE ENCLOSURE TO AVOID BUILD-UP OF EXPLOSIVE GAS CONCENTRATIONS**
- **ALWAYS** make an independent safety assessment of the experimental setup - Recommendations includes
  - Always to follow your own guide lines for safety
  - Always to follow ATEX rules/guidelines
  - Assessment of the amount of gas produced
  - Risks for reaching explosive gas limits
  - Always assume that all fittings/bottles etc are not gas tight
  - Measures to prevent/remove personal injuries if worst case scenario happens

Version 3.0 – Dec 22, 2024

# Flow rates & operation principle




After the cell, the liquid with H<sub>2</sub>/O<sub>2</sub> gases enters into the gas separation unit at the top and exits at the bottom, where after the liquid from the two separators enters into a common electrolyte bottle (typical for alkaline electrolysis with porous separators). The flow rate from the gas separation unit into the bottle is determined by (1) the liquid height difference between the bottle and the height in the gas separation unit (2) length and hydraulic resistance of the tubes between the gas separation unit and the bottle (3) Viscosity/temperature of the liquid and (4) the 'suction of the pump', head space of the bottle and how leak tight all fittings are. As an example lower left graph shows the liquid height inside the separation unit vs the flow rate for 1.6 mm ID tubes and 20 cm length (between separation unit and bottle) as function of flow rate. In general a relatively high liquid height is wanted as it ensures efficient gas separation, on the other hand too high levels will lead to an over-flow in the unit.



To ensure that the unit can work in many different configurations (flow rate, temperature etc) it comes with tubes (intended as connection between gas separation unit and bottle) that has both 1.6 mm ID and 2.2 mm ID. Going from 1.6 mm ID to 2.2 mm has a non-linear (exponential) increase of the hydraulic conductivity, while the length of the tube has a linear dependence. I.e. use of this unit requires some experimental work to calibrate the flow rate.

## Overview of included components

- 1 x Gas separation unit in PP (with all fittings and ferrules)
- 2 x 60 cm PP pipes for gas separation (OD  $\varnothing$ 16 mm ID  $\varnothing$ 13 mm) - for longer heights please make an inquiry
- 1 x Y fitting for 1/8 OD tubing in PEEK (including fittings/ferrules etc) - enables full integration with 'Hydraulic Connection Package for Flow Cells (Fittings, Lids, Tubings, Bottles) 1/8" '
- 2 m PTFE tube (1/8" OD & 1.6 mm ID)
- 2m PTFE tube (1/8" OD & 2.2 mm ID)




## Assembly







## Details on separation unit



H<sub>2</sub> exit (vent to a safe place)



## Mounting of tube inside pipe

